

Event 0:

- do nothing

Event 1:

- take picture
 - Turn on LED
- Aftermath: nothing

Event 2 (Solution pump 1 sec):

- Turn solution pumps on
 - Inject dextran + dye solution from one pump
 - Inject PEG solution from the other pump
- Wait for 1000 ms (PAUSE)
- Solution pumps off
- GOSUB Event 1

Event 3 (Pump reverse 1 sec):

- Turn solution pumps reverse
 - Retract dextran + dye solution from one pump
 - Retract PEG solution from the other pump
- Wait for 1000 ms (PAUSE)
- Solution pumps off
- GOSUB Event 1

Event 4 (Acid pump 1 sec):

- Turn acid pump on
- Wait for 1000 ms (PAUSE)
- Acid pump off

Event 5: Loop everything

- Set to Day 0 and Hour 29
- Replay events 1-4 (7 one second injection intervals followed by 7 one second retraction intervals)
- Write to event table

Loop 7 one second interval procedure

- Reset to right before event n occurs
 - Write to DurDay, DurHour, DurMin
 - Use TempReg to write to event table
 - \$FF code
 - **Alternatively manually 7 times using event table (refer to last year code)**
 - Boolean variable idea - times change before events, toggle variables

Event 6: Loop Event 1 36 times every 10 minutes

- Use counter variable to count times jumped back to DurDay 0 DurHour 23 DurMin 0
- If $i < 36$ then set (DurDay,DurHour,DurMin) = (0, 23, 59) and $i++$
- If $i \geq 36$ do nothing

1. Start experiment by waiting for a day
 - a. DurDay : 0
 - b. DurHour: 0
 - c. DurMin: 0
 - d. Event # : 0
2. Take Photo #1
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 0
 - d. Event # : 1
3. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 3
 - d. Event # : 2
4. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 6
 - d. Event # : 2
5. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 9
 - d. Event # : 2
6. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 12
 - d. Event # : 2
7. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 15
 - d. Event # : 2
8. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 18
 - d. Event # : 2
9. Begin pumping solution
 - a. DurDay : 1
 - b. DurHour: 0

- c. DurMin : 21
 - d. Event # : 2
- 10. Take Photo
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 24
 - d. Event # : 1
- 11. Loop Event 1
 - a. DurDay: 1
 - b. DurHour: 0
 - c. DurMin: 33
 - d. Event #: 6
- 12. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 36
 - d. Event # : 3
- 13. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 39
 - d. Event # : 3
- 14. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 42
 - d. Event # : 3
- 15. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 45
 - d. Event # : 3
- 16. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 48
 - d. Event # : 3
- 17. Reverse pumping
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin : 51
 - d. Event # : 3
- 18. Reverse pumping
 - a. DurDay : 1

- b. DurHour: 0
 - c. DurMin : 54
 - d. Event # : 3
- 19. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 0
 - c. DurMin: 57
 - d. Event # : 4
- 20. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 0
 - d. Event # : 4
- 21. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 3
 - d. Event # : 4
- 22. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 6
 - d. Event # : 4
- 23. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 9
 - d. Event # : 4
- 24. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 12
 - d. Event # : 4
- 25. Acid Pumps
 - a. DurDay : 1
 - b. DurHour: 1
 - c. DurMin: 15
 - d. Event # : 4
- 26. Loop everything
 - a. DurDay: 1
 - b. DurHour: 1
 - c. DurMin: 18
 - d. Set DurDay: 0
 - e. Set DurHour: 23

- f. Set DurMin: 0
- g. Event #: 5